

STEREOMICROSCOPE

MULTI-FUNCTION



Product Usage



| | |
|--------------------------|--|
| Weight | 4.9kg |
| Objective lens | Continuous zoom 7x~45x |
| Eyepiece | Binocular wide angle WF10X |
| Power | 220v-50hz, 110v-50/60hz |
| Transformer | Built-in 0-12v, continuous dimming |
| Light source on the lamp | 12v/15w halogen lamp |
| Under light source | 12v/15w halogen lamp |
| Eyepiece holder | 45° inclination, 360° horizontal rotation 55-75mm pitch |
| Packing dimensions | 39.2 X 24.8 X 37cm |

▼ CHARACTERISTICS

The microscope is suitable for a wide range of scenarios. Use a fully coated eyepiece to increase light transmission. A pure copper objective ensures image stability. Adjustable binoculars for different distances.

| Eyepiece | Standard configuration | | Additional large objective | | | | | |
|----------|------------------------|---------------|----------------------------|---------------|------------------|---------------|------------------|---------------|
| | | | 0.5X | | 1.5X | | 2X | |
| | Working distance | 100mm | Working distance | 165mm | Working distance | 45mm | Working distance | 30mm |
| | Magnification | Field of view | Magnification | Field of view | Magnification | Field of view | Magnification | Field of view |
| 10X/20 | 7X | 28.6 | 3.5X | 57.1 | 10.5X | 19 | 14X | 14.3 |
| | 45X | 4.4 | 22.5X | 8.9 | 67.5X | 3 | 90X | 2.2 |



Fully coated eyepiece
for increased light transmission

Suitable for various scenarios



LED halogen lamp




Variable power handle

Pure copper objective to ensure image stability





ORIGINAL PHOTO





Adjustable binocular tube

SHIPPING LIST



1.EYEPIECE PARTS

3.10/20X EYEPIECE X2

5.THE RACK

7.CHASSIS

9.THICKENED LARGE FOAM BOX

2. FOCUSING BRACKET

4. EYEPIECE COVER

6.INSTRUCTIONS

8.THICK LARGE CARTON